

Shropshire Continence Prescribing Guide – quick reference grid for catheter related items

SUPPLIER	PRODUCT	ORDER CODE	SIZE	PRODUCT	ORDER CODE
CATHETERS				CATHETERS	
Teleflex (Rusch)	Sympacath aquaflate female (hydrogel coated latex)	DH210112	12ch	Sympacath aquaflate male (hydrogel coated latex)	DH310112
	12 week licence	DH210114	14ch		DH310114
		DH210116	16ch	12 week licence	DH310116
		DH210118	18ch		DH310118
Teleflex (Rusch)	Brillant Aquaflate female (silicone catheter)	DA210112	12ch	Brillant Aquaflate male (silicone catheter)	DA310112
	12 week licence	DA210114	14ch		DA310114
		DA210116	16ch	12 week licence	DA310116
		DA210118	18ch		DA310118
Teleflex (Rusch)	PTFE Aquaflate female	DP210112	12ch	PTFE Aquaflate male	DP310112
	28 day licence	DP210114	14ch	28 day licence	DP310114
		DP210116	16ch		DP310116
		DP210118	18ch		DP310118
LEG BAGS				LEG BAGS	
Prosys	Long length – T Tap	P350L	350ml	Short length –T Tap	P350S
		P500L	500ml		P500S
		P750L	750ml		P750S
	Long Length - Lever Tap	P500L-LT	500ml	Short Length – Lever Tap	P500S-LT

Printed June 2018

Shropshire Continence Prescribing Guide – quick reference grid for catheter related items

SUPPLIER	PRODUCT	ORDER CODE	SIZE	PRODUCT	ORDER CODE
CATHETERS				CATHETERS	
Teleflex (Rusch)	Sympacath aquaflate female (hydrogel coated latex)	DH210112	12ch	Sympacath aquaflate male (hydrogel coated latex)	DH310112
	12 week licence	DH210114	14ch		DH310114
		DH210116	16ch	12 week licence	DH310116
		DH210118	18ch		DH310118
Teleflex (Rusch)	Brillant Aquaflate female (silicone catheter)	DA210112	12ch	Brillant Aquaflate male (silicone catheter)	DA310112
	12 week licence	DA210114	14ch		DA310114
		DA210116	16ch	12 week licence	DA310116
		DA210118	18ch		DA310118
Teleflex (Rusch)	PTFE Aquaflate female	DP210112	12ch	PTFE Aquaflate male	DP310112
	28 day licence	DP210114	14ch	28 day licence	DP310114
		DP210116	16ch		DP310116
		DP210118	18ch		DP310118
LEG BAGS				LEG BAGS	
Prosys	Long length – T Tap	P350L	350ml	Short length –T Tap	P350S
		P500L	500ml		P500S
		P750L	750ml		P750S
	Long Length - Lever Tap	P500L-LT	500ml	Short Length – Lever Tap	P500S-LT

Printed June 2018

Shropshire Continence Prescribing Guide – quick reference grid for catheter related items

SUPPLIER	PRODUCT	ORDER CODE	SIZE	PRODUCT	ORDER CODE
	NIGHT BAGS			NIGHT BAGS	
Prosys	2 litre drainable (single use)	PSU2	2000ml	2 litre drainable (7 day use if continually attached)	P2000
	SECURING STRAP			SECURING STRAP	
Coloplast	G Strap – standard (50cm)	383001	-----	G strap – abdominal 90cm	383003
	CATHETER VALVE			PROSYS LEG BAG SLEEVE	
Coloplast	Simpla Catheter Valve	CV3808	-----	Small 24-40cms (yellow)	PLS3881
Bard	Flip- Flo catheter valve	BFF5	-----	Medium 36-50cms (blue)	PLS3904
			-----	Large 40-65cms (red)	PLS3928
	ANAESTHETIC GEL			ANAESTHETIC GEL	
Optimum medical	Optilube 6ml For urethral & suprapubic	1160	6ml or 11 ml	Optilube 11ml For urethral & suprapubic	1161
Teleflex	Cathejell Mono 8.5g for pts allergic to lidocaine	CJM8501	8.5g or 12.5g	Cathejell mono 12.5g for pts allergic to lidocaine	CJM12501
	CATHETER MAINTENANCE				
Bbraun	Urotainer NaCL (0.9%)	FB99849	50ml	To check patency, remove blood clots or debris	
	Urotainer Suby G (3.23%)	9746609	2x30ml	For internal encrustation alongside pH monitoring	
	Urotainer Solution R (6%)	9746625	2x30ml	For external catheter encrustation	

Printed June 2018

Shropshire Continence Prescribing Guide – quick reference grid for catheter related items

SUPPLIER	PRODUCT	ORDER CODE	SIZE	PRODUCT	ORDER CODE
	NIGHT BAGS			NIGHT BAGS	
Prosys	2 litre drainable (single use)	PSU2	2000ml	2 litre drainable (7 day use if continually attached)	P2000
	SECURING STRAP			SECURING STRAP	
Coloplast	G Strap – standard (50cm)	383001	-----	G strap – abdominal 90cm	383003
	CATHETER VALVE			PROSYS LEG BAG SLEEVE	
Coloplast	Simpla Catheter Valve	CV3808	-----	Small 24-40cms (yellow)	PLS3881
Bard	Flip- Flo catheter valve	BFF5	-----	Medium 36-50cms (blue)	PLS3904
			-----	Large 40-65cms (red)	PLS3928
	ANAESTHETIC GEL			ANAESTHETIC GEL	
Optimum medical	Optilube Active 6ml For urethral & suprapubic	1160	6ml or 11 ml	Optilube Active 11ml For urethral & suprapubic	1161
Teleflex	Cathejell Mono 8.5g for pts allergic to lidocaine	CJM8501	8.5g or 12.5g	Cathejell mono 12.5g for pts allergic to lidocaine	CJM12501
	CATHETER MAINTENANCE				
Bbraun	Urotainer NaCL (0.9%)	FB99849	50ml	To check patency, remove blood clots or debris	
	Urotainer Suby G (3.23%)	9746609	2x30ml	For internal encrustation alongside pH monitoring	
	Urotainer Solution R (6%)	9746625	2x30ml	For external catheter encrustation	

Printed June 2018